

Contents

<i>Preface</i>	<i>v</i>
<i>Contributors</i>	<i>ix</i>
1 Setting Up an Arabidopsis LED Culture Module that Simulates Plant Neighbor Proximity	1
<i>Carlos D. Crocco</i>	
2 Photobody Detection Using Immunofluorescence and Super-Resolution Imaging in Arabidopsis	7
<i>Giorgio Perrella, Anna Zioutopoulou, Andrew Hamilton, and Eirini Kaiserli</i>	
3 Analysis of Shade-Induced Hypocotyl Elongation in Arabidopsis	21
<i>Yetkin Çaka Ince and Vinicius Costa Galvão</i>	
4 Isolation of UVR8 Protein Complexes	33
<i>Yan Liu and Xi Huang</i>	
5 Phenotypic Study of Photomorphogenesis in Arabidopsis Seedlings	41
<i>Chuanwei Yang, Famin Xie, and Lin Li</i>	
6 Experimental Procedures for Studying Skotomorphogenesis in <i>Arabidopsis thaliana</i>	49
<i>Huanhuan Jin, Hong Li, and Ziqiang Zhu</i>	
7 Global Identification for Targets of Circadian Transcription Factors in Arabidopsis and Rice Using Chromatin Immunoprecipitation Followed by Sequencing (ChIP-seq)	61
<i>Shuxuan Xu, Jing Huang, Jian Jin, and Wei Huang</i>	
8 Co-immunoprecipitation Assays to Detect In Vivo Association of Phytochromes with Their Interacting Partners	75
<i>Pengyu Song, Shaoman Zhang, and Jigang Li</i>	
9 Detection of UVR8 Homodimers and Monomers by Immunoblotting Analysis in Tomato	83
<i>Guoqian Tang, Xiaorui Liu, and Li Lin</i>	
10 Characterization of Seedling Greening Process in Plant Photomorphogenesis	95
<i>Wanqing Wang, Yuhong Li, and Rongcheng Lin</i>	
11 Protoplast System for Studying Blue-Light-Dependent Formation of Cryptochrome Photobody	105
<i>Xiangguang Lyu, Hongyu Li, and Bin Liu</i>	
12 Uncover the Nuclear Proteomic Landscape with Enriched Nuclei Followed by Label-Free Quantitative Mass Spectrometry	115
<i>Tan Wang, Zhuang Lu, and Lei Wang</i>	

Contents

13	Strategies to Study Dark Growth Deficient or Slower Mutants in <i>Chlamydomonas reinhardtii</i>	125
	<i>Huanling Yang, Fei Han, Yue Wang, Wenqiang Yang, and Wenfeng Tu</i>	
14	Co-immunoprecipitation Assay for Blue Light-Dependent Protein Interactions in Plants.	141
	<i>Jingyi Zhang and Shengbo He</i>	
15	Detecting Blue Light-Dependent Protein-Protein Interactions by LexA-Based Yeast Two-Hybrid Assay	147
	<i>Xiaolong Hao and Ling Li</i>	
16	Express Arabidopsis Cryptochrome in <i>Sf9</i> Insect Cells Using the Baculovirus Expression System	155
	<i>Xu Li, Tawen Liu, and Hongtao Liu</i>	
17	Semi-In-Vivo Pull-Down Assay for Blue Light-Dependent Protein Interactions	161
	<i>Xu Li, Tawen Liu, and Hongtao Liu</i>	
18	Tobacco System for Studying Protein Colocalization and Interactions.	167
	<i>Jingyi Zhang and Shengbo He</i>	
	<i>Index</i>	175